

CytoCHECK SPAchip[®] pH Green Single-Detection Kit

CytoCHECK SPAchip[®] pH Green Single-Detection Kit allows measurement of **intracellular and extracellular pH levels** by changes in green fluorescence intensity, which allows a more comprehensive study of the living single-cell physiology and provides additional insights using imaging analyzers.

SPAchip[®] technology allows long-term intracellular monitoring without cytotoxic effects, ensuring minimal interference with cellular processes. Designed for living single-cell analysis, this kit provides researchers with a powerful, ready-to-use tool to study key biological processes, including metabolism, cell viability, and response to environmental changes, while maximizing experimental reproducibility and efficiency.

Highlights

- Non-toxic for living single cells. Measures intracellular and extracellular pH levels by changes in fluorescence intensity.
- Live and quantitative, long-term monitoring of intracellular and extracellular pH changes.
- Composed of fluorescently labeled silicon microparticles that can be internalized in the cytosol of cultured cells.
- Provides a more comprehensive study of single cell physiology and metabolism.
- Compatible with most imaging platforms.
- Cell type flexibility.
- Ready-to-use, robust workflow.




Figure 1: SH-SY5Y cell line (neuroblastoma cells) fluorescence images containing internalized SPAchip[®] pH Green for intracellular and extracellular pH detection. Nuclei stained in blue and actin cytoskeleton in red are visualized alongside SPAchip[®] pH-sensitive green fluorescence. Images A, B, and C represent different fields of view, showcasing the distribution and internalization of SPAchips[®] within living cells. This technology enables live and quantitative monitoring of intracellular pH, making it an essential tool for studies in cell metabolism, drug development, and toxicity assessment.

Contact us
a4cell.com | info@a4cell.com

Arrays for Cell Nanodevices SL

CAIT - UPM. Campus de Montegancedo s/n, 28223 Pozuelo de Alarcón – Madrid (Spain)
+34 910 679 519 - info@a4cell.com

Product features

- SPAchip® assay kits are novel cell-based assays for living single-cell that bring together the fields of **nanotechnology and cell biology**.
- CytoCHECK SPAchip® kits are composed of **fluorescently labeled silicon microparticles - SPAchips®**- that can be internalized in cultured cells to monitor changes in specific intracellular analyte concentrations for long periods of time.
- **CytoCHECK SPAchip® pH greenSingle-Detection Kit** allows measurement of **intracellular and extracellular pH levels** by changes in fluorescence intensity.
- **High Sensitivity and Stability:** unlike traditional pH-sensitive dyes, CytoCHECK SPAchip® pH green Single-Detection Kit ensures stable, reproducible pH measurements over extended periods of time, making it ideal for applications in phenotypic screening, drug discovery, and toxicity assessment.
- **Optimized for Imaging and Flow Cytometry:** compatible with fluorescence microscopy, high-content screening (HCS), and flow cytometry, enhancing imaging-based analysis for intracellular and extracellular pH studies.



Each CytoCHECK SPAchip pH Green Single-Detection Kit contains:

~2.5x10⁶
ASSAY SPAchips

ASSAY SPAchip® tube (embedded in a solid fluorescence-protective soluble film)

5 mL

ASSAY buffer tube (Sterile, cell culture suitable)

~2.5x10⁵ CONTROL SPAchips/100 µL

CONTROL SPAchip® tube (non-fluorescent, ready-to-use)

CytoCHECK SPAchip® pH Green Single-Detection Kit

Product code	S-001-PHG
Amount	~2.5 millions of SPAchips
Applications	Cell viability, proliferation, cell image acquisition
Assay time	30 minutes
Assay type	Living single-cell based
Analyte	pH
Detection method	Green fluorescence*
Fluorescence	λ _{ex} : 488 nm; λ _{em} : 520 nm
Measuring range	pH 4.5 - 7.5
Compatible Platforms	Fluorescence microscopy, HCS/HCA platforms (recommended 20x magnification and over) and flow cytometry
Sample type	Adherent cells, suspension cells
Solubility	Soluble in assay buffer (aqueous)

*Ensure to follow the full User Protocol

Contact us

a4cell.com | info@a4cell.com

Arrays for Cell Nanodevices SL

CAIT - UPM. Campus de Montegancedo s/n, 28223 Pozuelo de Alarcón – Madrid (Spain)

+34 910 679 519- info@a4cell.com

